

| | |
|--|----------|
| DamID-seq protocol | 1 |
| Tissue preparation | 1 |
| DNA extraction (Zymo Quick-DNA Miniprep plus kit) | 1 |
| DpnI digestion | 1 |
| Adaptor ligation | 2 |
| DpnII digestion | 2 |
| PCR amplification | 2 |
| QC and sonication | 3 |
| Library preparation | 3 |
| Purify with SeraMag beads | 3 |
| Adjust concentrations | 4 |
| End repair | 4 |
| Adenylate 3' Ends | 4 |
| Adaptor ligation | 4 |
| Clean-up with SeraMag beads | 4 |
| Clean-up with SeraMag beads again | 5 |
| Enrich DNA fragments | 5 |
| Clean-up with SeraMag beads (final time) | 5 |
| Library QC and multiplexing | 6 |
| Bioinformatic processing of samples | 6 |
| Notes | 6 |
| Citations | 6 |
| Reagents and equipment | 7 |
| DamID buffers | 8 |
| Homebrew sequencing buffers | 8 |
| Ligation of DamID adaptors | 9 |
| DamID oligos | 9 |
| Sequencing adaptors and primer sequences | 10 |
| Annealing adaptors | 10 |
| Notes re. adaptors | 10 |
| Preparation of SeraMag beads | 11 |
| Test beads by purifying Bioline Hyperladder V | 11 |
| Notes | 12 |
| Optional: homemade magnetic rack for bead purification | 12 |
| Troubleshooting | 13 |
| Secondary peak after PCR | 13 |
| Using the tapestation | 14 |
| Using the Qbit | 14 |

DamID-seq protocol

This protocol is based on the following published research:

Marshall OJ, Southall TD, Cheetham SW, Brand AH (2016) *Nat Protocols* 11:1586–1598

Marshall OJ, Brand AH (2015) *Bioinformatics* 31:3371-3

Southall TD, Gold KS, Egger B, Davidson CM, Caygill EE, Marshall OJ, Brand AH (2013) *Dev Cell* 26:101-12

Vogel MJ, Peric-Hupkes D, van Steensel B (2007) *Nat Protoc* 2:1467–1478.

Tissue preparation

Aim to process enough material to give >10,000 driven cells

1. Take embryos, dissected tissue or heads from -80°C freezer (stored without buffer)
2. Add 75 µL **sterile H₂O** + 20µL **500mM EDTA**
NB: Addition of EDTA is optional, but recommended when working with tissue containing the gut (or any tissue with high concentrations of nucleases/proteases). We have not tried lower concentrations than 100mM; lower concentrations may also work.
3. **If working with embryos or heads**, use sterile (wash in 100% ethanol) pestle to homogenise samples (smush heads by rotating pestle between fingers and pressing down to bottom of tube)

DNA extraction (Zymo Quick-DNA Miniprep plus kit) (2hrs)

[Day 1]

The following assumes the use of a vacuum manifold. If using without a manifold, add in centrifugation steps as per the manufacturer's protocol.

1. Mix together 10 µL **Proteinase K** (Zymo Quick-DNA kit) + 95µL **Solid Tissue Buffer** (Zymo Quick-DNA kit; blue) per sample and vortex to make a **Prot K master mix**.
2. Add 105µL of the **Prot K master mix** to each sample and flick gently to mix; digest at 56°C for 1-3 hours (heat block)
3. Cool to RT, add 400µL **Genomic Binding Buffer** (Zymo Quick-DNA kit) and mix by gentle inversion
4. Add all solution to a spin column on vacuum manifold and draw through
5. Add 400µL **DNA Pre-Wash Buffer** (Zymo kit) and draw through
6. Add 700µL **g-DNA Wash Buffer** (Zymo kit) and draw through
7. Add 200µL **g-DNA Wash Buffer** and draw through
8. Transfer column to a collection tube and spin at maximum speed in benchtop centrifuge for 2 mins to dry the column
9. Transfer spin column to a new 1.5 ml tube, add 50 µL **DNA Elution Buffer** and leave at RT for at least 30 min. Spin at >6000 x g for 1 min and keep the flow-through
10. Optional: run 1µl on a 1% gel to check quality (should be a single band on the top of gel and not a smear)

DpnI digestion (overnight + 30 mins)

1. Transfer 43.5µl of elution to a new 1.5ml tube
2. Add 5µL CutSmart Buffer and 1.5µl DpnI, mix very gently and digest at 37°C overnight.
(This digestion step can be reduced to two hours if required.)

DpnI can be optionally heat-inactivated for 20 mins at 80°C.

[Day 2]

3. Clean-up with Machery-Nagel PCR Purification kit (100µL **Buffer NTI**; 2x 700µL washes); elute in 32µL DEPC H₂O [**Label tube "A"**]

4. Transfer 15µL of sample into an 8-well PCR strip; store the remainder at -20°C as a backup. *There is generally very little/undetectable amounts of DNA at this point, as uncut genomic DNA (which should be the majority of DNA in most use cases) will not pass through the spin columns – the purified DNA should only come from induced cells.*
*(If working with large amounts of tissue and/or if concerned about DNA shearing, the concentration of samples can be measured and a maximum of 750ng should be used going forward. However, such a high yield is highly unusual. **Very low or undetectable yields at this stage produce excellent data – do not throw samples away!**)*

Adaptor ligation (2.5 hrs)

1. Add 4µL **Adaptor Ligation Buffer (AL)**
2. Add 1µL **T4 DNA ligase** and mix well
3. Ligate for 2 hours (16°C for 2 hours → 65°C for 10 mins in PCR machine)

DpnII digestion (2.5 hrs)

1. Add 19µL **DpnII Digestion Buffer (DIIB)**
2. Add 1µL **DpnII enzyme** and mix well
3. Digest at 37°C for 2 hours
4. Heat inactivate at 65°C for 20 mins

PCR amplification (2.5 hrs + 30 mins)

1. Add 118µL **PCR Buffer (PCR)**
2. Add 2µL **MyTaq polymerase** and mix well
3. Split into 4 x 40µL reactions in a PCR strip using multichannel
4. Run PCR (use DamID 17 programme in PCR machine):
 - 68°C 10 min
 - 94°C 30 sec
 - 65°C 5 min
 - 68°C 15 min
 - Repeat 3x
 - 94°C 30 sec
 - 65°C 1 min
 - 68°C 10 min
 - Repeat 17x
 - 94°C 30 sec
 - 65°C 1 min
 - 68°C 2 min
 - 68°C 5 min
 - 4°C hold

This PCR protocol is from Vogel et al., 2007. Both the initial long extension for 10 mins and the subsequent long cycles appear to be required. A 3 min extension time, as present in the final 17 cycles, should amplify all fragments below 5kb; these represent 99.97% of all GATC fragments in the D. melanogaster genome.

(An argument could be made that fragments larger than 5kb will not represent genuine binding at that fragment, given Dam-fusion proteins typically only methylate within +/- 2kb of their binding site in a seemingly normal distribution. On that basis, longer amplification times should probably be discouraged.)

5. PCR purify (160µL sample + 320µL **Buffer NTI**). [Day 3]
If using MyTag, wash column three times with kit's EtOH wash buffer before elution. Use 32µL DEPC H₂O to elute (leave water on column for at least 5 min) [**Label tubes "C"**]
6. Run 1µL on a gel to check quality (should be a smear between ~400bp and 2kb)
7. Spec DNA conc on nanodrop

QC and sonication (>30 mins + overnight)

1. Dilute samples: 2µg DNA in 90µL DEPC H₂O in 1.5mL TPX sonication tubes
Using dedicated sonication tubes is important to obtain consistent results
2. Add 10µL **CutSmart Buffer** and mix well; cool on ice
3. Sonicate in Diagenode Bioruptor Nextgen:
 - a. Add 500mL ice to chamber, top up to the line with water and mix
 - b. Cool tube holder cassette on ice before using
 - c. Sonicate for 5 cycles (27 secs on, 27 secs off) on high power
The aim is to achieve a ~300bp average fragment size. All sonicators are different; if using a different sonication device you will need to optimise the conditions first.
4. Check fragment size on Tapestation (genomic tape) – ensure that average size is ~300bp
Slight variations in fragment size are acceptable. Very large fragments (>600bp) may impede clustering efficiency and sequencing yields
5. Add 1µL **AlwI enzyme** and digest overnight at 37°C
AlwI removes the DamID adaptors and initial GATC sequence (cuts at GGATCNNNN|N). Digestion can be carried out either before or after sonication; it makes no difference to the end result. AlwI cannot be heat inactivated, but will not cause issues if carried over in subsequent steps.
6. Transfer 70µL (or 90µL) of each sample into PCR 8-well strips for library prep

Library preparation (~6 hrs)

[Day 4]

The following protocol uses homebrew reagents and is not compatible with Illumina kits. If you wish to use a commercial Illumina kit, we recommend the TruSeq Nano kit – follow the manufacturer's protocol except for the PCR step, in which case still use a reduced number of cycles (6 cycles seems optimum for most cases – see notes below). Commercial AmpureXP beads are a (vastly more expensive) alternative to homemade Seramag beads.

Purify with Seramag beads (40 mins)

1. Add 105µL of Seramag beads to 70µL sample (or 135µL to 90µL) and mix well
Pipette mix 20x and then flick mix the strip and pulse down. Rapid mixing of beads and sample is important. Use 90µL of samples if concentration prior to sonication was low.
2. Incubate at RT, 10'
3. Place on magnetic stand 10' (or until clear)
4. Remove supernatant
5. Wash twice in 190µL 80% EtOH (30 secs wash time)
6. Let plate stand for 5' to air dry
7. Resuspend in 25µL **Resuspension Buffer** and remove from magnetic stand
8. Mix well and incubate 2' RT
9. Place on magnetic stand for 5' (or until clear)
10. Put 22.5µL of the supernatant into a new, clean tube for the next step

Adjust concentrations (30 mins)

1. Measure library conc on Qubit
2. Dilute samples to (no more than) 500ng of DNA in 20 μ L **Resuspension Buffer**

End repair (1 hr)

1. Add 7.5 μ L **End Repair Buffer (ERB)**
2. Add 2.5 μ L **End Repair enzyme** and mix well
3. Run PCR program “Truseq 1 2 ER Hi”:
 - a. Incubate for 30 mins at 30°C
 - b. Heat inactivate enzymes for 20 mins at 75°C

Adenylate 3' Ends (40 mins)

1. Add 0.75 μ L **Klenow 3' to 5' exo-** enzyme and mix well
2. Run PCR program “Truseq 2 AT”:
 - a. Incubate for 30 mins at 37°C
3. Proceed immediately to adaptor ligation

Adaptor ligation (20 mins)

1. Add 2.5 μ L **NEB Quick Ligase** enzyme
2. Add 2.5 μ L of relevant adaptor
3. Run PCR program “Truseq 3 lig”:
 - a. Incubate for 10 mins at 30°C
4. Add 5 μ L **Stop Ligation Buffer**

Clean-up with Seramag beads (35 mins)

1. Add 40 μ L Seramag beads and mix well
2. Incubate at RT, 10'
3. Place on magnetic stand 5' (or until clear)
4. Remove supernatant
5. Wash twice in 190 μ L 80% EtOH (30 secs wash time)
6. Let plate stand for 5' to air dry
7. Resuspend in 52.5 μ L **Resuspension Buffer** and remove from magnetic stand
8. Mix well and incubate 2' RT
9. Place on magnetic stand for 5' (or until clear)
10. Put 50 μ L of the supernatant into a new, clean tube for the next step

Clean-up with Seramag beads again (35 mins)

1. Add 50µL Seramag beads and mix well
2. incubate at RT, 10'
3. Place on magnetic stand 5' (or until clear)
4. Remove supernatant
5. Wash twice in 190µL 80% EtOH (30 secs wash time)
6. Let plate stand for 5' to air dry
7. Resuspend in 22.5µL **Resuspension Buffer** and remove from magnetic stand
8. Mix well and incubate 2' RT
9. Place on magnetic stand for 5' (or until clear)
10. Put 20µL of the supernatant into a new, clean tube for the next step

This second clean-up is required to ensure complete removal of sequencing adaptor dimers – should these be present in the final library they will vastly out-compete genuine library sequences when hybridising to the flow cell.

Enrich DNA fragments (20 mins)

1. Add 5µL **PCR Primer Cocktail**
2. Add 25µL **PCR Master Mix**
3. Run PCR program “Truseq 4 PCR 6 cycles”:
98°C for 30 secs
6 cycles of:
 98°C for 10 secs
 60°C for 30 secs
 72°C for 30 secs
72°C for 5 mins

Note: aim to avoid concatamer formation (see troubleshooting section below)

Clean-up with Seramag beads (final time) (35 mins)

1. Add 50µL Seramag beads and mix well
2. incubate at RT, 10'
3. Place on magnetic stand 5' (or until clear)
4. Remove supernatant
5. Wash twice in 190µL 80% EtOH (30 secs wash time)
6. Let plate stand for 5' to air dry
7. Resuspend in 32.5µL **Resuspension Buffer** and remove from magnetic stand
8. Mix well and incubate 2' RT
9. Place on magnetic stand for 5' (or until clear)
10. Put 30µL of the supernatant into a new clean tube

Library QC and multiplexing (1 hr)

1. QC sample with tapestation (use genomic tape) – ensure that concatemers are not present (see “Troubleshooting” below) and that adapter dimers are not present (small peak around 240bp) in the sample. Record the average fragment length for each sample.
2. Measure DNA conc with Qubit (see attached protocol)
3. Calculate molarity of each sample from fragment length and concentration – use the spreadsheet template (ask Owen for this), or see the “Using the Qbit” section for details as to how to calculate molarity manually
4. Pool samples to give 20nM final DNA conc (with all libraries at equal concentration) – use the spreadsheet template to calculate this automatically
5. Sequence as single-end 50nt (SE50) on Illumina sequencer.

Bioinformatic processing of samples

1. $\text{Log}_2(\text{Dam-fusion/Dam})$ binding profiles can be obtained from either raw reads in FASTQ format or aligned BAM files using the `damidseq_pipeline` software, freely available from https://owenjm.github.io/damidseq_pipeline. Detailed usage instructions and examples are available from the website.
2. Profiles can be viewed using browser software such as IGV, as detailed on the website.
3. Calling expressed genes from RNA pol II ratio files generated by the `damidseq_pipeline` software can be carried out using the `polii.gene.call.r` script, freely available from <https://github.com/owenjm/polii.gene.call>.

Notes

1. Wear gloves, use filter tips and clean, autoclaved H_2O to avoid contamination
Exceptions to filter tips: DNeasy/Qiagen kits before final elution
Always use filter tips when preparing sequencing libraries
The use of DEPC H_2O throughout this protocol, while technically unnecessary, reflects the need to use a clean, contamination-free source of H_2O
2. **Be very gentle with samples prior to DpnI digestion.** Any genomic DNA shearing will result in broken ends that can potentially ligate DamID adapters. This will lead to random signals. Do not vortex, but mix samples gently by inversion and/or slow pipetting with a blue tip (and/or *very* gentle flicking). Always check the quality of the DNeasy-prepped DNA on a gel before continuing with the protocol.
3. For forward planning, each sample uses 278.5 μL beads during the sequencing library prep

Citations

If you find this protocol useful, please cite

<http://www.ncbi.nlm.nih.gov/pubmed/23792147> (original TaDa publication) and
<http://www.ncbi.nlm.nih.gov/pubmed/26112292> (bioinformatic processing and NGS method); or
<http://www.ncbi.nlm.nih.gov/pubmed/27490632> (Nature Protocols)

Reagents and equipment

MyTaq HS DNA polymerase (BioLine, Cat# BIO-21112 (200 μ l) or BIO-21113 (500 μ l))

Note: it's been hard to find a replacement for the old Clontech Advantage cDNA polymerase. MyTaq HS seems the best so far; Clontech's Advantage 2 will work, but at reduced yield and is not recommended. Another, as yet untested, alternative would be Cesium Klentaq AC LA from DNA Polymerase Technology.

Seramag SpeedBead, 3 EDAC/PA5, 1 μ m, 15mL (Fisher Scientific, Cat# 12326433)

Note: Seramag beads need to be diluted in PEG before use as described in the dedicated section below. AmpureXP beads (which are ready to use and do not need further preparation) may be substituted for Seramag beads.

NEBNext® High-Fidelity 2X PCR Master Mix 1.25mL (NEB, Cat# M0541S)

*Alternative to NEBNext PCR MM (works equally well, but may be more expensive):
Kapa HiFi HS ReadyMix 1.25mL (Anachem, Cat# KK2601)*

Quick Ligase kit (NEB, Cat# M2200S)

We're just using the Quick ligase enzyme, which is simply NEB T4 ligase at 2000 U/ μ L. However, at \$97 per 30 μ L, this is cheaper than high conc. NEB T4 ligase (#M0202T) at \$64 per 10 μ L.

NEB enzymes:

T4 ligase (NEB, Cat# M0202S)
T4 DNA polymerase (NEB, Cat# M0203S)
T4 polynucleotide kinase (NEB, Cat# M0201S)
Klenow fragment (NEB, Cat# M0210S)
Klenow 3' to 5' exo- (NEB, Cat# M0212L)

NEB restriction enzymes:

DpnI (Cat# R0176L)
DpnII (Cat# R0543L)
AlwI (Cat# R0513L)

dNTPs (NEB, Cat# N0447S)

Kits

Zymo Quick-DNA miniprep plus (Cat# D4069)

Macherey-Nagel NucleoSpin Gel and PCR Clean-up (Cat# 740609.250)

Other equivalent kits also work; we find these to be the most effective in terms of price, performance and ease of use. Note that for genomic DNA extraction, we recommend a kit that allows efficient elution in 50 μ L.

DamID buffers

| Ligation buffer (AL) | | x100 | Aliquot |
|--------------------------------------|------|-------------|----------------|
| 10x T4 ligase buffer | 2 | 200 | 40 |
| 0.8uL ds ADR | 0.8 | 80 | |
| H2O | 1.2 | 120 | |
| DpnII digestion buffer (DIIB) | | x100 | Aliquot |
| Dpn II Buffer | 4 | 400 | 190 |
| H2O | 15 | 1500 | |
| PCR | | x50 | Aliquot |
| 5X MyTaq HS Buffer | 32 | 1600 | 1180 |
| DamID-PCR primer (50 μ M) | 2.5 | 125 | |
| DEPC H2O | 83.5 | 4175 | |

Homebrew sequencing buffers

The following buffers and components are mostly adapted from [Ethan Ford's protocol](#)

| End Repair Buffer (ERB) | 1x (7.5μL) | Master mix x50 |
|--------------------------------|----------------------------------|-----------------------|
| NEB T4 ligase buffer 10x | 3.0 | 150 |
| 10mM dNTPs | 1.2 | 60 |
| DEPC H ₂ O | 3.3 | 165 |
| Total | 7.5 | 375 |

| End Repair enzymes | 1x (2.5μL) | Master mix x50 |
|---|----------------------------------|-----------------------|
| T4 DNA polymerase (NEB, 3U/ μ L) | 1.14 | 56.82 |
| Klenow Fragment (NEB, 5 U/ μ L) | 0.23 | 11.36 |
| T4 polynucleotide kinase (NEB, 10 U/ μ L) | 1.14 | 56.82 |
| Total | 2.5 | 125.00 |

PCR Master Mix (PCR MM)

NEBNext 2x HiFi Hotstart Ready-mix (aliquot into 105 μ L aliquots and use as supplied)

PCR Primer Cocktail (PPC)

25 μ L PCR1 primer (100 μ M)
 25 μ L PCR2 primer (100 μ M)
 50 μ L DEPC H₂O

Resuspension Buffer (RSB)

10mM Tris-HCl, pH 8.0
 0.1 mM EDTA

Ligation Stop Solution

0.5M EDTA

Ligation of DamID adaptors

Making ds AdR stock (50uM):

1. Take 50 μ L AdRt (100uM) and 50 μ L AdRb (100 uM)
2. Incubate in removable metal heating block, 95°C, 2 mins
3. Remove heating block and allow to cool to room temperature (should take > 45mins)

DamID oligos

AdRt

CTAATACGACTCACTATAGGGCAGCGTGGTCGCGGCCGAGGA

AdRb

TCCTCGGCCG

Annealed structure:

```
5' CTAATACGACTCACTATAGGGCAGCGTGGTCGCGGCCGAGGA 3'
                                     ||| ||| ||| |||
3' GCCGGCTCCT 5'
```

DamID_PCR

GGTCGCGGCCGAGGATC

Sequencing adaptors and primer sequences

(adaptor barcodes highlighted in red; additional bases for higher indices in blue)

| | |
|------------------|--|
| Universal | AATGATACGGCGACCACCGAGATCTTACACTCTTTCCCTACACGACGCTCTTCCGATC*T |
| Index 1 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATCAGG ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 2 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGATGT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 3 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTAGGC ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 4 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TGACCA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 5 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ACAGTG ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 6 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCCAAT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 7 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CAGATC ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 8 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ACTTGA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 9 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GATCAG ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 10 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TAGCTT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 11 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GGCTAC ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 12 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTTGTA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 13 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGTCAACA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 14 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGTTCCGT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 15 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATGTCAGA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 16 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CCGTCCG ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 18 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GTCCGCAC ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 19 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GTGAAACG ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 20 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GTGGCCTT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 21 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GTTTCGGA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 22 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGTACGTA ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 23 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GAGTGGAT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 25 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ACTGATAT ATCTCGTATGCCGTCTTCTGCTT*G |
| Index 27 | /5Phos/GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATTCCTTT ATCTCGTATGCCGTCTTCTGCTT*G |
| PCR1 | AATGATACGGCGACCACCGA*G |
| PCR2 | CAAGCAGAAGACGGCATAACGA*G |

/5Phos/ = 5' Phosphorylation; * = Phosphorothioate linkages

We order adaptors from IDT @ 100 nmol synthesis + standard desalting

Oligonucleotide sequences © 2015 Illumina, Inc. All rights reserved. Derivative works created by Illumina customers are authorized for use with Illumina instruments and products only. All other uses are strictly prohibited.

Annealing adaptors

1. Resuspend adaptor oligos at 100µM in TE + 50mM NaCl
2. Mix 25µL of relevant index + 25µL Universal primer in PCR tube
3. Incubate in removable metal heating block, 95°C, 2 mins
4. Remove heating block and allow to cool to room temperature (should take > 45mins)

Notes re. adaptors

1. If multiplexing four or fewer libraries, selecting adaptors with barcodes that are too similar may result in a reduced number of reads passing the filter. In this case, the preferred indexes to use (in order) are 4, 7, 6 and 8.


Preparation of Seramag beads

(modified from Rohland N and Reich D. (2012) Cost-effective, high-throughput DNA sequencing libraries for multiplexed target capture. *Genome Res* 22:939–46. doi:10.1101/gr.128124.111)

1. In 15mL Falcon tube add:
 - 3g PEG-8000 (20% w/v final conc)
 - 3mL 5M NaCl (1M final)
 - 150 μ L 1M Tris pH 8.0 (10mM final)
 - 30 μ L 0.5M EDTA (1mM final)
2. add DEPC H₂O to 14mL
3. mix by inversion until PEG is dissolved (~5mins)
4. mix Speedbeads container well to ensure beads are in suspension
5. add 300 μ L Seramag Speedbeads and mix by inversion
6. add DEPC H₂O to 15mL volume and store at 4°C

Test beads by purifying Bioline Hyperladder V

1. Make mastermix of diluted ladder: 10 μ L ladder + 20 μ L DEPC H₂O per sample
2. Purify ladder with 1x, 1.25x, 1.5x amounts of beads
3. Elute in 20 μ L DEPC H₂O
4. Run on tapestation using D1K tape (should look like 20% PEG-8000 lanes):



Notes

1. Rohland and Reich used 18% PEG-8000; however, this does not appear to perform as well as 20% in ladder purification (compared with AmpureXP)
2. Based on current prices:
 Ampure beads are £689.67 for 60mL working solution
 Seramag beads are £317.62 for 15mL, which makes 750mL working solution
 Seramag beads are therefore more than 27x more cost-effective (!)


Optional: homemade magnetic rack for bead purification

Thermoscientific provide a 96-well magnetic plate for £763, which works reasonably well. However, the total retail cost of the neodymium magnets in this plate is ~£6 (!)

Whilst it would be possible to buy individual ring magnets and glue them to a 96-well rack with epoxy resin, I've found it easier to make a simple 8-well (i.e. one PCR strip) magnetic plate using:

- 1x N42 Neodymium magnet (75x10x3mm) (e.g. [these magnets](#))
- 1x 20 μ L filter tip box insert (acting as the 96-well support)
- 1x 75x10x2.5mm plastic spacer (supplied with a pair of the strip magnets)

The magnet together with spacer is taped to the tip box insert as shown below.



Important: as the magnet needs to remain firmly in contact with the edge of the strip tubes at all times, use a second PCR strip and/or tape on the side of the strip to hold the magnet in place.


The advantage of using a strip magnet rather than a circular magnet is that the beads are pulled to one side of the tube rather than sitting in a ring, greatly facilitating liquid removal. (The disadvantage is that such a configuration is not compatible with 96-well high-throughput solutions).

Strength N42 seems to work well, but there should be no harm in using more powerful magnets if you can obtain them (the only result will be a faster purification process).

Troubleshooting

Secondary peak after PCR

This generally only happens using 8 or 10 (or more) cycles as PCR reaction is exhausted – a secondary peak twice the size of the original peak is seen on the tapestation genomic plots. This peak results from concatemers of the amplified product. If you see this, reduce number of PCR cycles and try again (6 cycles seems to work well for TaDa).



Using the tapestation

1. Prepare ladder and samples in 8-well PCR strips:
 - 2µL pre-made ladder dilution in lane 1
 - 3µL sample buffer + 1µL sample in subsequent lanes
 - (Pre-made ladder dilution: ladder:loading-dye 1:10 or 1:15. Using 1:15 will save significant reagent costs but will obviously affect quantitation.)*
2. Cap strip and mix well by vortexing (pulse down and mix several times)
3. Turn on tapestation and laptop
4. Launch Tapestation software and wait until prompted to insert the tape
5. Open tapestation, place tape in slot in correct orientation (barcode towards you on the right)
 - We recommend using genomic tapes rather than DIK HS tapes – the genomic tapes give a better picture of adaptor-concatemer formation if present. D5K tapes may present an even better solution here, but we haven't tried these yet.*
6. Insert strip with ladder and samples and close the machine
7. Select the wells you want to run on the software, and add labels in the table that appears
8. Click “Start” to begin the run (will take ~10 min to complete)
9. Once run is complete and analysis software is open, use the region settings to determine average fragment sizes and save a report

Using the Qbit

1. Prepare samples as per the included protocol.
2. NB – use 1µL sample + 199µL working solution for each sample
3. With each sample, press the button to obtain the DNA concentration
4. DNA molarity $\cong \frac{1,500}{\text{fragment size in bp}} \times [\text{concentration in ng}/\mu\text{L}]$
5. Use a spreadsheet template to calculate final pooling of samples

Qubit™ Assays

www.invitrogen.com/qubit

QUICK
REFERENCE
CARD

NOTE: For best results, store the dye and the buffer at room temperature. Store the DNA, RNA, and protein standards at 4°C. Ensure that all assay reagents are at **room temperature** before you begin.

1. Set up two Assay Tubes for the standards (three for the protein assay) and one tube for each user sample.
2. Prepare the Qubit™ **Working Solution** by diluting the Qubit™ reagent 1:200 in Qubit™ buffer. Prepare 200 µL of **Working Solution** for each standard and sample.
3. Prepare the Assay Tubes* according to the table below.

| | Standard Assay Tubes | User Sample Assay Tubes |
|--|----------------------|-------------------------|
| Volume of Working Solution (from step 2) to add | 190 µL | 180–199 µL |
| Volume of Standard (from kit) to add | 10 µL | — |
| Volume of User Sample to add | — | 1–20 µL |
| Total Volume in each Assay Tube | 200 µL | 200 µL |

* Use only thin-wall, clear 0.5 mL PCR tubes. Acceptable tubes include Qubit® assay tubes (set of 500, Cat. no. Q32856) or Axygen PCR-05-C tubes (VWR, part no. 10011-830).


4. Vortex all tubes for 2–3 seconds.
5. Incubate the tubes for 2 minutes at room temperature (15 minutes for the Qubit™ protein assay).
6. Insert the tubes in the Qubit® 2.0 Fluorometer and take readings. For detailed instructions, refer to the Qubit® 2.0 Fluorometer manual.
7. *Optional:* Using the Dilution Calculator feature of the Qubit® 2.0 Fluorometer, determine the stock concentration of your original sample.

— MP39808 | MAN0002474 —
Rev. date: 28 Sept 2010

molecular
probes® | invitrogen™
by life technologies™

For technical support, email tech_support@invitrogen.com.

For country-specific contact information, visit www.invitrogen.com.



* where n = number of Standards plus number of Samples

For research use only. Not intended for any animal or human therapeutic or diagnostic use.
©2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.



Vortex all assay tubes for 2-3 seconds

Incubate at room temperature for 2 minutes (15 minutes for Qubit™ protein assay)

Read tubes in Qubit® 2.0 Fluorometer

